

DOSE RECORDS DATA MANAGEMENT SOFTWARE

Dose records data management software is designed to read **total dose** from the dosimeter automatically once the dosimeter is placed on the docking station. It also facilitates Configuration of Alarm Settings for Cumulative Dose and Dose-rate, assigning Unique ID for individual dosimeters, Calibration of the dosimeters for Cumulative Dose. The Application software routes data from PC to Dosimeter through virtual comm. port connecting to Docking station. This in turn routes the data through IR port to the Dosimeter. Once the data is retrieved into PC, it gets recorded in a database. Users can retrieve and analyze dose records data personal ID-wise & date wise.

Docking station with this software is an ideal choice for institutions who may like to maintain dose records of the personnel working in a Radiological installation or facility.

